

Seldinger Guide Wire



Dimensions

Ø 0.018" / 0,46 mm
Ø 0.021" / 0,53 mm
Ø 0.025" / 0,63 mm
Ø 0.028" / 0,71 mm
Ø 0.032" / 0,81 mm
Ø 0.035" / 0,89 mm
Ø 0.038" / 0,96 mm

Length

30 – 450 cm

Surface

Stainless steel
PTFE coated

Distal tip form

Straight tip
J-tip radius: 1,5 – 3 – 6 – 15 mm

Distal flexible length

Standard 4,5 cm flexible tip


Shaft stiffness

Standard
Heavy duty
Customised

Mode of delivery

Sterile – CE marked
Bulk OEM

Customised products on request



SP Medical develops, manufactures and provides plastics and surface solutions to the global medical devices and pharmaceutical industries. Our core capabilities are injection molding, surface treatment and processing of plastic products. At our state-of-the-art facilities in Denmark and Poland, we manufacture components, medical devices, single-use devices, and tablets with active pharmaceutical ingredients, and we apply function-enhancing coatings to a number of products.

SP Medical is a company in SP Group A/S, an international plastics group headquartered in Denmark. SP Group is listed on NASDAQ OMX Copenhagen.

www.sp-medical.com